

NEOPULSER DESIGN

The Neopulser is based on the work of Dr. Raymond Rife and Dr. Hulda Clark, pioneers in electro-energetic medicine. This design produces proprietary square-wave frequencies from 15Hz – 150,000 KHz which will produce upper harmonics that are clustered closer together and thus hit more frequencies in the lower range. These etheric tools include: an [orgone generator](#) composed of a plastic resin/metal chip matrix material sometimes called **orgonite**, a **Mobius coil**, **Amethyst & Garnet** crystals, a strong **neodymium magnet**, and a **containment coil**.

The Neopulser puts out a fast-rise, positive-offset, **pulsing DC** waveform, which is extremely rich in harmonics. The sudden shock of this burst waveform can cause **scalar waves** to be formed in the presence of a Mobius coil. Scalar waves can transfer information of a complex nature and can transfer the 'frequency signature' of orgone, the crystals, etc., into the body. The orgone generator increases the effectiveness of the etheric or subtle energy level. Most people simply don't know that the frequency signatures of both healthy and unhealthy body components are written or encoded within the etheric body (the astral body), which surrounds and permeates the physical body. The energetic signature of disease conditions first shows up in the etheric body and will then translate into physical manifestation if the condition is not addressed in sufficient time. This reality has been established and proven repeatedly by Kirlian photography, Radionics, the work of Cleve Backster, etc.

DIRECTIONS:

Open the back of the unit, insert a 9-Volt battery and close up the unit. When we ship your order the battery is already included and placed in the device. Click the plastic toggle over and it will blink green. This indicates that the unit is ON. The green light must be on for the unit to work. All units are tested at our facility before going out. If the green light does not come on, replace the battery. If the device still doesn't work, please email us: Products@ascendedhealth.com

Note: When we ship your order the battery is already included and placed in the device

Our recommendation is to use the device on the sole of the foot in a sock while resting or sleeping or in the palm of the hand with or without a glove. The copper discs must be touching bare skin to benefit from the pulses. Leave the unit on for as long as you feel comfortable. **IF YOU PLACE THE UNIT ANYWHERE ELSE ON THE BODY MAKE SURE TO MOVE IT AROUND. IF YOU DO NOT MOVE THE UNIT, YOU MAY GET A SLOW BURN ON YOUR SKIN.** Unless you use a pad with salt water (see instructions below) or if the copper is touching the sole of the foot or palm of the hand.

GETTING STARTED:

It's a good idea to begin your pulsing experience by wearing the Neopulser on the **sole of your foot** when you go to bed at night. This will give you six to eight hours of "zapping time" without irritation, reddening, or darkening skin. You may also build up your time more each day, adding 30 more minutes daily. If you are in good health, you won't get much of a die-off reaction and you probably won't notice any symptoms. Make sure to keep your body hydrated by drinking proper amount of water while you are using this device. If you have a more serious health condition, you need to balance the amount of time you zap and the effort you take to detoxify and clear the die-off waste. Drink water often, eat light foods (especially freshly prepared veggies, fruit, & juicing), and take detoxifying baths (such as Mustard Bath). This device works synergistically with the Ascended Health Anti-Venom detoxifying formula and Triple Skin Oil as well as ProAlive and Active Detox probiotics.

We suggest that you pulse around the clock for the first month. It is the constant application of electricity to the body that seems to have a negative effect on the various parasites, venom, bacteria, and viruses in the body. Our body can take it, but the simple organisms cannot. During sleep, it works very well – this is the time your body takes to repair itself.

DETOXIFYING EFFECTS:

Do not pulse for any LESS THAN 20 minutes at first. This may stop the complete clearing of the whole size-range of parasites and could result in sickness. The effectiveness of this unit is based on the low-grade electricity you can send through your body for *long periods at a time*. Our human bodies can take on electricity, but parasites, bacteria, and viruses (in our opinion) cannot -- they need to escape or die.

We have never heard of anyone being harmed using the Neopulser or similar devices. Some people may experience momentary dizziness, nausea, gas, or headache from the parasitic die-off. Worms, for instance, have bacteria, viruses, and fungi within them. They are harmless to the worm but could be harmful to the human if they are not cleared as soon as they are released from the worm that's no longer alive.

The Neopulser also seems to drive stored toxins, including toxic metals, out of the tissues for elimination. This is done, we believe, through the stimulation of Kundalini (55hz) frequencies, which is another form of cellular energy.

FIELD EFFECTS:

The Neopulser also produces a vibrational energy field effect and can speed up healing and pain relief when worn close to the affected area. The field effect of the Neopulser extends about a foot in all directions.

SKIN REACTIONS:

Some people may get some sort of skin reaction after wearing the Neopulser for a while. The copper from the discs can react with the salt from your sweat and produce **copper chloride**, which can irritate the skin. The skin reddening will eventually stop as the blood turns more alkaline. In this case, our guess is that the body is compensating in some way and produces a buffering compound to neutralize the irritant. **We have not seen this happen to the soles of the feet.**

If you cover the discs with 100% cotton pads or cotton flannel or cotton batting saturated with [a 10-25% solution of saltwater](#) can eliminate the skin eruptions completely. The saltwater is now carrying the current to the skin without direct contact with the copper, avoiding the electrolysis problem. The only downside is that you have to keep the cotton damp so it can conduct electricity and carry the signal into the body.

Using the device on the soles of the feet at night will prevent any skin eruption even without using the cotton pads. The palms of the hands will also resist skin eruptions. The skin cells on the soles of the feet and the palms of the hands are of a special nature that will not react to copper metal electrolysis.

NOTE: Ascended Health's 60 days refund policy does not apply to the Neopulser device. Thank you.